
Faichney Center Termómetro Clínico Digital Dual Scale

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
Conforme a ASTM E 1112

Faichney Medical Company
17517, Fabrica Way, Suite G
Cerritos, CA 90703, USA

1122307 - Rev 04
Faichney Dual Scale Thermometer
Reorder No. 6210 / 6211 / 6212 / 6213

Faichney Dual Scale Digital Clinical Thermometer

PLEASE READ CAREFULLY BEFORE USING

BODY TEMPERATURE

Normal body temperature can vary from 97°F (36.1 °C) to 99 °F (37.2 °C) . The generally accepted "Normal" temperature is 98.6 °F (37.0 °C). Usually your temperature is lower in the morning than in the afternoon. Normal temperature varies according to the location of the temperature reading. While the oral temperature follows the guidelines previously described, a rectal temperature is generally 1.0 °F (0.5 °C) higher. Conversely, an axillary (under the arm) temperature will be 1.0 °F (0.5 °C) lower.

GENERAL CARE AND USE

BEFORE USE PLEASE NOTE

1. The thermometer should be held during use. The patient should be seated or lying down.
2. TEMPERATURES SHOULD BE TAKEN UNDER ADULT SUPERVISION.

SPECIAL PRECAUTIONS

3. Avoid exposure to temperature extremes, dust, excessive humidity and direct sunlight.
4. Do not soak the thermometer in water or other liquids. Keep parts other than the probe away from moisture.
5. Do not attempt to sterilize with heat or by boiling.
6. Do not bend or drop the unit. It is not shock proof.
7. Do not allow the unit to come into contact with any thinners or chemical solvents.
8. Do not attempt to disassemble the unit except to replace the battery.
9. Do not bite the main body or sensor unit.

DIRECTIONS FOR CLEANING

Clean the Thermometer probe before and after use by wiping it with soap and water. Rinse with 70% Isopropyl Alcohol, then rinse thoroughly in cold water.

SPECIFICATIONS

DISPLAY RESOLUTION	1/10 DEGREE FOR BOTH °C & °F
MEASUREMENT RANGE	89.6 °F to 109.4 °F (32.0 °C to 43.0 °C)
ACCURACY	±0.3 °F / ±0.2 °C (89.6 °F to < 95.0 °F / 32 °C to < 35 °C) ±0.2 °F / ±0.1 °C (95.0 °F to < 104.0 °F / 35 °C to 40 °C) ±0.3 °F / ±0.2 °C (>104.0 °F to 109.4 °F / >40 °C to 43 °C)
MEASUREMENT METHOD	PEAK HOLD SYSTEM
BATTERY	ONE 1.5V DC ALKALINE MANGANESE BUTTON SIZE LR41 TYPE OR EQUIVALENT
BATTERY LIFETIME	UPTO 3 YEARS UNDER NORMAL USE AND STORAGE TEMPERATURE OR (27 °F)
AUTOMATIC SHUT-OFF	AFTER APPROXIMATELY 5 MINUTES
OPERATING TEMPERATURE	16 °C to 40 °C (61°F to 104 °F)
RELATIVE HUMIDITY	15% to 95% NON CONDENSING
STORAGE TEMPERATURE	-20 °C to 50 °C (- 4 °F to 122 °F)
RELATIVE HUMIDITY	15% to 95% NON CONDENSING

CAUTION






NEVER SUBMERGE THIS THERMOMETER IN ANY LIQUID.
NEVER ATTEMPT TO STERILIZE WITH HEAT, BY BOILING OR IN A GAS OR STEAM AUTOCLAVE.

Thermometer performance may degrade in the following cases :

- (a) Operation/storage outside the specified operating/storage temperature range respectively.
- (b) Mechanical Shock/dropping.
- (c) Patient temperature is below ambient temperature.

OPERATING INSTRUCTIONS

1. Ensure that the thermometer has been turned off before taking temperature.
2. The unit is factory set to °F. To start the unit in °F, Press the blue colored ON/OFF switch momentarily once. To start the unit in °C, Press and hold the blue colored ON/OFF switch for approximately 2 sec., when °F start blinking press the switch again one to choose °C. If you need to go back of °F scale, press and hold ON/OFF switch for approximately 2 sec and choose your scale by momentary press of switch, once °C starts blinking.
3. The following appear in quick sequence.
 - (a) This is the LCD display test (See Fig. 1). It appears for approximately 2 seconds and indicates all the segments of the LCD are functioning properly.
 - (b) If the previous reading is between 89.6 °F ~ 109.4 °F, (32°C ~ 43°C), the unit will display the last time measured temperature value for about 3 seconds with symbol M (Fig.2.). If the last time measured temperature value is out of range (< 32 °C & > 43 °C OR < 89.6 °F & > 109.4 °F) the unit will display 97.7 °F (36.5 °C) without the M symbol (Fig.3) for about 3 seconds.
 - (c) The unit goes through the internal self checking test procedure. If the unit is good, it will display 98.6 °F (37.0 °C) (Fig.4) for about one second, otherwise it will display "Err" (Fig.5)
 - (d) When "Lo" appears with the flashing °F (°C) sign (Fig. 6), the unit is ready to take temperature.

ORAL USE

1. Close mouth for about 5 minutes before measurement.
2. Place the thermometer tip under the tongue on a heat pocket (✓) as shown in Fig. 7.
3. Hold the thermometer in the same spot for about 40 seconds.
4. Readings will continue to increase and the °F (or °C) symbol will flash during measurement (See Fig.8).
5. When peak temperature is reached, the °F (or °C) symbol will stop flashing (See Fig. 9) and
 - a) When the temperature is higher than 99.4 °F (or 37.4 °C) the beeper will sound "beep-beep-beep" ten times to indicate a high temperature.
 - b) When the temperature is lower than or equal to 99.4 °F (or 37.4 °C) the beeper will sound "beep" ten times.
6. After measurement, record the temperature, then turn off the unit to save the battery life. If you forget to turn it off, the unit will power down automatically after approximately 5 minutes.

Fig.7

Fig. 8

Fig. 9

NOTE : THE MOUTH MUST REMAIN CLOSED DURING THE MEASUREMENT TO ENSURE ACCURATE READINGS.

RECTAL USE

Rectal temperature should only be taken when oral and axillary temperatures are inconvenient. It is recommended that axillary temperatures be taken on infants. LUBRICATE WITH LUBRICATING JELLY AND INSERT THE PROBE NO MORE THAN 1/2 INCH INTO RECTUM. Hold thermometer in place until you hear the beep signal and the °F (or °C) sign stops flashing (about 40 seconds.)


AXILLARY USE

Dry armpit before measurement. Place the tip of the probe under arm so that the tip is touching the skin. This area is called axilla. Hold the thermometer so that it is parallel to the patient's body and fold the arm over the chest. This ensures that the tip of the probe is well covered and not influenced by air. Hold thermometer in place until you hear the beep signal and the °F (or °C) sign stops flashing (about 3 minutes)

BATTERY REPLACEMENT

Replace battery with an alkaline LR 41 or equivalent type whenever the low battery indicator "±" flashing.


1. Remove screw located on the back of the thermometer. Press down on the battery cover and push in the direction show by the arrow in the illustration.
2. Remove battery with a non-metallic probe.
3. Insert new battery. NOTE : The +Vé side must face match with "+" sign in plastics.
4. Replace the cover.
5. Dispose off used battery properly.

NOTE : Inserting battery with reverse polarity may damage the unit. Dispose off this thermometer properly as per local regulations. This Digital clinical thermometer is pre-calibrated and can not be re-calibrated.

DISPOSABLE PROBE COVER USE.

1. Insert the thermometer into probe cover. Take care not to push the thermometer tip through the end of the probe cover.
2. Holding the thermometer and tab of the probe cover as shown, pull away the covering longer side.
3. Tilt thermometer to remove cover from shorter side.
4. Dispose off the probe cover after taking temperature.

NOTE : For rectal use, after applying probe cover lubricate thermometer tip with a water soluble jelly.

LIFETIME GUARANTEE

Your new Faichney Dual Scale Digital Clinical Thermometer is guaranteed to be free from manufacturer's defects for life under normal use.

However, if for any reason (other than abuse or normal battery replacement) your thermometer does not function properly we will either repair or replace it free of charge for life. This warranty is valid only if unit returned for repair is properly packed to avoid damage in shipping. All parts are covered by the warranty, except batteries and probe covers, where applicable.

Please check the battery before returning it. This is usually the source of the problem, should repair be needed, ship freight paid to :

Conforms to ASTM E 1112

Faichney Medical Company
17517, Fabrica Way, Suite G
Cerritos, CA 90703, USA

1122307 - Rev 04
Faichney Dual Scale Thermometer
Reorder No. 6210 / 6211 / 6212 / 6213